



**IMPLANTS DIFFUSION
INTERNATIONAL**

Developer and Manufacturer
since 1987

WWW.IDI-DENTAL.COM



ID'BOX

ID-SPASSOR

The solution to manage the volume of the prosthetic tooth



Step 1 :

Use the drill adapted to the diameter of the future crown, stick it to the side of the adjacent tooth and proceed to drilling.



Step 2 :

Remove the drill and replace by the parallelism guide of the same color.



The **ID-SPASSOR** guide includes :

- 6 drilling guides (the drills are 7 mm length and Ø2mm)
- 6 assorted non-drilling guides (virtual teeth) with diameter from 4,5mm to 12mm.

