


# IMPLUS MC



**NOTE:** The letter "S" identifies implants with STRAIGHT NECK. For these implants, countersink is NOT necessary. Only for implants with diameter 3,75 we recommend the use of the countersink in case of bone quality D1 - D2.

3.6

Unique Platform



	Implant Ø3,75 mm	Implant Ø4,1 mm	Implant Ø4,5 mm	Implant Ø5,5 mm
<b>Core Ø</b>	3,0 mm	3,3 mm	3,5 mm	4,5 mm
<b>Apex Ø</b>	2,3 mm	2,5 mm	2,7 mm	3,7 mm
<b>Thread pitch</b>	0,6 mm	0,6 mm	0,75 mm	0,75 mm
<b>Neck Ø</b>	4 mm	4,1 mm	4,5 mm	5,5 mm
<b>Implant length (h)</b>	<b>Code</b>	<b>Code</b>	<b>Code</b>	<b>Code</b>
<b>8,0 mm</b>	01IMC3708TTS	01IMCS4108TTS	01IMCS4508TTS	01IMCS5508TTS
<b>10,0 mm</b>	01IMC3710TTS	01IMCS4110TTS	01IMCS4510TTS	01IMCS5510TTS
<b>11,5 mm</b>	01IMC3711TTS	01IMCS4111TTS	01IMCS4511TTS	01IMCS5511TTS
<b>13,0 mm</b>	01IMC3713TTS	01IMCS4113TTS	01IMCS4513TTS	01IMCS5513TTS
<b>15,0 mm</b>	01IMC3715TTS	01IMCS4115TTS	01IMCS4515TTS	//